

IN-VITRO BIOEQUIVALENCE

How do you prove a drug works the same when the body never absorbs it?

Sucroferic Oxyhydroxide 500 mg chewable tablets · Phosphate binder



A study-design walkthrough

Swipe →

THE PROBLEM

No absorption means **no blood levels to compare**

1

A binder, not a systemic drug

The active, polynuclear iron(III)-oxyhydroxide, is practically insoluble. It works inside the gut and leaves the body unabsorbed.

2

Classic PK studies don't apply

With nothing reaching the bloodstream, a standard pharmacokinetic BE study can't show sameness.


→

So the test moves to the lab

USFDA's product-specific guidance asks generics to show equivalence through **in-vitro phosphate binding studies**.

HOW IT WORKS

It traps dietary phosphate by ligand exchange



Binding has to hold across the gut's whole pH range, which is why the in-vitro study tests at two very different conditions:

pH 3.0

Stomach-like, acidic

pH 7.5

Intestine-like, near neutral

WHAT USFDA ASKS FOR

Two binding studies, run under tight control

E Equilibrium

At least 8 phosphate concentrations, binding measured at equilibrium

K Kinetic

At least 8 time points up to 24 hours, showing how fast it binds

- **Whole tablets** incubated, at pH 3.0 and pH 7.5
- **37 °C** with constant gentle shaking
- **pH checked every 15 minutes** and adjusted until equilibrium holds
- **12 replicates** of every binding study

THE SCALE OF IT

One "simple" binding claim. Thousands of data points.

8

pivotal studies: 2 equilibrium +
6 kinetic (low, mid, high at
each pH)

96

sets: every study repeated 12
times

3,840

samples in the pivotal phase
alone

~7,000

injections across validation
(~2,500) and pivotal (~4,500)

16

parameters per method
validation, to ICH Q2 and
USFDA

~3,650

test and reference tablets, end
to end

FROM THE BENCH

Two tablets of the same drug don't behave the same in a flask

1 Dispersion

One product dispersed in 2–3 minutes; the other needed 6–8 minutes with sonication.

2 pH drift

Starting pH and the acid needed to hold it differed widely, before and after phosphate spiking.

3 Phase separation

Settling times ranged from 5 to 15 minutes, which changes when you can sample.



So we optimised every variable

Sonication time, shaker RPM and pH control, with ~100% recovery confirmed in aqueous and placebo controls.

INSIDE THE LAB

The study, up close




1 Whole tablets dispersing



2 pH adjustment, before phosphate



3 Final samples on the shaker



4 Phase separation, test vs reference

TIMELINE

From first experiment to submission-ready data

2

MONTHS

Method development

Two methods, one per binding study

2

MONTHS

Method validation

Two validations, at pH 3.0 and pH 7.5

1

MONTH

Pilot studies (optional)

Eight studies to de-risk the pivotal run

3

MONTHS

Pivotal studies

Eight studies × 12 replicates, test vs reference

Total, end to end

~8 months



DOCTORS'
ANALYTICAL LABORATORIES PVT. LTD.
Testing for a better world

Planning an in-vitro BE package for a **non-absorbed drug?**

Our pharma division in Navi Mumbai does:

In-Vitro Binding Studies

NG/G Tube Studies

IVRT / IVPT

Method Development & Validation

Multipoint Dissolution



US FDA Inspected



Health Canada Approved